

For - BEVERAGE & FOOD / COSMETIC / PHARMACEUTICALS / CHEMISTRY / METROLOGY

## PBDT-1 PET bottle drop tester

PBDT-1 PET bottle drop tester, allows a repeatable drop testing. Referenced to the ISBT requirements, bottles up to 160 mm in diameter, can be dropped from 800 to 1950 mm height AT2E PBDT-1 includes also an angle adjustment for special test.

The walls are movable and adjustable to fit any kind of bottles from 150 ml to 4 L  $\,$ 

PBDT-1 allows testing in vertical or horizontal position.

Safe conception and stainless steel structure ensure durability and easy cleaning.

Draining output is installed under the bottom plate.

A front door gives an easy access to see the test and retrieve the bottles after impact.

## Advantages:

- Vertical or Horizontal drop test
- · Test repeatable with height adjustment ruller
- Transparent safety door for easy observation and removal sample
- · Visible test and fall
- stainless steel design safe and easy to clean
- Position limit baffle adjustable, compatible with products with different sizes

## **Specifications:**

Bottle range : round or square from 150ml to 4L

Max bottle diameter: 160mm

Max bottle height: 400mm

Max Drop Height: 1950 mm

Min Drop Height: 800 mm

Dimensions: 950mm x 960mm x 2050mm (W x L x H)

Weight: 250 kg

## **Optional Parts:**

Model with angle adjustment





info@at2e.com Tél.: +33 (0)1 64 66 03 02